

NEW SIZES!
NOW INDIVIDUALLY PACKAGED

Richmond

**Soft
BiteBlock™**

AVAILABLE IN:
Large
Medium
Small



- Reduces the potential of compression on the endotracheal tube or laryngeal mask airway
- Place between the back molars to protect the teeth & tongue - pressure is exerted on the back molars instead of the front teeth
- Individually packaged to protect the product from contamination until time of use (bulk is not individually packaged)
- 100% medical-grade cotton
- Latex free & absorbent
- Saves labor, cost and time (versus manually rolling gauze)



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
260000	SOFT BITEBLOCK™	3/4" DIA. X 4", BULK, LARGE, NON-STERILE	120/BOX	3 BOXES	USA
260004	SOFT BITEBLOCK™	3/4" DIA. X 4", LARGE PACKAGED, NON-STERILE	100/BOX	3 BOXES	USA
260003	SOFT BITEBLOCK™	3/4" DIA. X 2.75", MEDIUM PACKAGED, NON-STERILE	100/BOX	3 BOXES	USA
260005	SOFT BITEBLOCK™	9/16" DIA. X 2.5", SMALL PACKAGED, NON-STERILE	100/BOX	3 BOXES	USA

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.

Richmond
RichmondDentalAndMedical.net
855-841-0743



Manufactured where grown, in the USA for Richmond, Charlotte, NC

7/20/18