

# Endotracheal Tube

Artificial Ventilation Channel




TUORen



## Meet your Endotracheal Tube Needs

Tuoren is a healthcare solutions provider focusing on medical consumables and extending to electronic medical devices, surgery devices and biomedical materials. International subsidiaries are established in the United States and India. Its products are exported to over 60 countries and regions in the world.

The disposable endotracheal tube is suitable for inserting through the nose and mouth. It is suitable for clinical anesthesia, emergency patients and other patients who need mechanical ventilation.




# Endotracheal Tube

## ENDOTRACHEAL TUBE MADE OF PVC

Size: 2.0#-10.0#

- Available with cuff and without cuff.
- Cuffed tube with high volume and low pressure, less compression for airway.



### Cuff

Providing even pressure to maintain good sealing, reducing pressure on the tissues of the trachea

### Without cuff



CUFF			
SIZE	CAT. NO	SIZE	CAT. NO
2.0	ET0020	6.5	ET0065
2.5	ET0025	7.0	ET0070
3.0	ET0030	7.5	ET0075
3.5	ET0035	8.0	ET0080
4.0	ET0040	8.5	ET0085
4.5	ET0045	9.0	ET0090
5.0	ET0050	9.5	ET0095
5.5	ET0055	10.0	ET0100
6.0	ET0060	10.5	-


UNCUFF			
SIZE	CAT. NO	SIZE	CAT. NO
2.0	ET1020	6.5	ET1065
2.5	ET1025	7.0	ET1070
3.0	ET1030	7.5	ET1075
3.5	ET1035	8.0	ET1080
4.0	ET1040	8.5	ET1085
4.5	ET1045	9.0	ET1090
5.0	ET1050	9.5	ET1095
5.5	ET1055	10.0	ET1100
6.0	ET1060	10.5	-

# Reinforced Endotracheal Tube

*TUBE BODY MORE FLEXIBLE,  
EASY TO BE RESHAPED*

Size: 2.5#-10.0#

- Built-in stainless steel spring with the feature of anti-kinking.
- Black mark easy to operate.



**Reinforced tube body**

Built-in stainless steel spring with the feature of anti-kinking

CUFF			
SIZE	CAT. NO	SIZE	CAT. NO
2.0	-	6.5	RET0065
2.5	RET0025	7.0	RET0070
3.0	RET0030	7.5	RET0075
3.5	RET0035	8.0	RET0080
4.0	RET0040	8.5	RET0085
4.5	RET0045	9.0	RET0090
5.0	RET0050	9.5	RET0095
5.5	RET0055	10.0	RET0100
6.0	RET0060	10.5	-

UNCUFF			
SIZE	CAT. NO	SIZE	CAT. NO
2.0	-	6.5	RET1065
2.5	RET1025	7.0	RET1070
3.0	RET1030	7.5	RET1075
3.5	RET1035	8.0	RET1080
4.0	RET1040	8.5	RET1085
4.5	RET1045	9.0	RET1090
5.0	RET1050	9.5	RET1095
5.5	RET1055	10.0	RET1100
6.0	RET1060	10.5	-

Black mark easy to operate

# TPU-Cuffed Reinforced Endotracheal Tube

THE CUFF IS TRANSPARENT WITH BETTER ADHESION

Size: 2.5#-8.0#

- The cuff is more transparent and thinner.
- Better adhesion to the wall, easy to insert, especially for children.



Reinforced			
SIZE	CAT. NO	SIZE	CAT. NO
2.0	-	6.5	RETST0065
2.5	RETST0025	7.0	RETST0070
3.0	RETST0030	7.5	RETST0075
3.5	RETST0035	8.0	RETST0080
4.0	RETST0040	8.5	-
4.5	RETST0045	9.0	-
5.0	RETST0050	9.5	-
5.5	RETST0055	10.0	-
6.0	RETST0060	10.5	-

# Oral/Nasal Endotracheal Tube

## TUBE BODY IS PREMOLDED


Size: 2.0#-9.0#

- Suitable for oral, maxillofacial, nasal cavity and eye surgery.
- The tube is far away from the operation area, so as to have a better view of the operation.



Oral Endotracheal tube CUFF					
SIZE	CAT. NO	CAT. NO (REINFORCED)	SIZE	CAT. NO	CAT. NO (REINFORCED)
2.0	-	-	6.5	OET0065	OET0065
2.5	OET0025	OET0025	7.0	OET0070	OET0070
3.0	OET0030	OET0030	7.5	OET0075	OET0075
3.5	OET0035	OET0035	8.0	OET0080	OET0080
4.0	OET0040	OET0040	8.5	OET0085	OET0085
4.5	OET0045	OET0045	9.0	OET0090	OET0090
5.0	OET0050	OET0050	9.5	-	-
5.5	OET0055	OET0055	10.0	-	-
6.0	OET0060	OET0060	10.5	-	-

Oral Endotracheal tube UNCUFF					
SIZE	CAT. NO	CAT. NO (REINFORCED)	SIZE	CAT. NO	CAT. NO (REINFORCED)
2.0	-	-	6.5	OET1065	OET1065
2.5	OET1025	OET1025	7.0	OET1070	OET1070
3.0	OET1030	OET1030	7.5	OET1075	OET1075
3.5	OET1035	OET1035	8.0	OET1080	OET1080
4.0	OET1040	OET1040	8.5	OET1085	OET1085
4.5	OET1045	OET1045	9.0	OET1090	OET1090
5.0	OET1050	OET1050	9.5	-	-
5.5	OET1055	OET1055	10.0	-	-
6.0	OET1060	OET1060	10.5	-	-



Nasal Endotracheal tube CUFF					
SIZE	CAT. NO	CAT. NO (REINFORCED)	SIZE	CAT. NO	CAT. NO (REINFORCED)
2.0	-	-	6.5	NET0065	NET0065
2.5	NET0025	NET0025	7.0	NET0070	NET0070
3.0	NET0030	NET0030	7.5	NET0075	NET0075
3.5	NET0035	NET0035	8.0	NET0080	NET0080
4.0	NET0040	NET0040	8.5	NET0085	NET0085
4.5	NET0045	NET0045	9.0	NET0090	NET0090
5.0	NET0050	NET0050	9.5	-	-
5.5	NET0055	NET0055	10.0	-	-
6.0	NET0060	NET0060	10.5	-	-

Nasal Endotracheal tube UNCUFF					
SIZE	CAT. NO	CAT. NO (REINFORCED)	SIZE	CAT. NO	CAT. NO (REINFORCED)
2.0	-	-	6.5	NET1065	NET1065
2.5	NET1025	NET1025	7.0	NET1070	NET1070
3.0	NET1030	NET1030	7.5	NET1075	NET1075
3.5	NET1035	NET1035	8.0	NET1080	NET1080
4.0	NET1040	NET1040	8.5	NET1085	NET1085
4.5	NET1045	NET1045	9.0	NET1090	NET1090
5.0	NET1050	NET1050	9.5	-	-
5.5	NET1055	NET1055	10.0	-	-
6.0	NET1060	NET1060	10.5	-	-

# Injectable Endotracheal Tube

*HAVE INJECTION PORT TO INJECT DRUGS*

Size: 6.0#-8.0#

- Can inject different medicines according to the patient's condition.
- Reduce the occurrence of adverse patient reactions.

Injectable			
SIZE	CAT. NO	SIZE	CAT. NO
6.0	IET060	7.5	IET075
6.5	IET065	8.0	IET080
7.0	IET070		




# Supraglottic & Subglottic Injectable Endotracheal Tube

*TWO SPRAY CUFFS ARE SET ON THE UPPER AND LOWER GLOTTIC*

Size: 6.0#-8.0#

- The uniquely designed injection device allows the medicine to be sprayed onto the mucous membrane of the respiratory tract in a mist form through the hole. The medicine is widely distributed, safe and effective.

Injectable			
SIZE	CAT. NO	SIZE	CAT. NO
6.0	SIET060	7.5	SIET075
6.5	SIET065	8.0	SIET080
7.0	SIET070		





# Suction Endotracheal Tube

**HAVE SUNCTION CHANNEL  
MAINLY FOR THE SUCTION  
OF SUBGLOTTIC SECRETIONS**

Size: 6.0#-8.0#

- Have normal and reinforced type.
- Sputum suction incision reduces the occurrence of VAP.
- The operation is simple and convenient.



Normal			
SIZE	CAT. NO	SIZE	CAT. NO
6.0	SET0060	7.5	SET0075
6.5	SET0065	8.0	SET0080
7.0	SET0070	8.5	SET0085



Reinforced			
SIZE	CAT. NO	SIZE	CAT. NO
6.0	-	7.5	SRET0075
6.5	SRET0065	8.0	SRET0080
7.0	SRET0070	8.5	SRET0085

# BlockBuster Endotracheal Tube

## HIGH SUCCESS RATE FOR BLIND ENDOTRACHEAL INTUBATION

Size: 3.0#-8.0#

- Special soft flexible tip reduce tracheal wall trauma and subglottic mucosal lesions.
- The olive-shape cuff reduce the airway trauma.



BlockBuster Endotracheal Tube			
SIZE	CAT. NO	SIZE	CAT. NO
2.0	-	6.5	BRETE0065
2.5	-	7.0	BRETE0070
3.0	BRETE0030	7.5	BRETE0075
3.5	BRETE0035	8.0	BRETE0080
4.0	BRETE0040	8.5	-
4.5	BRETE0045	9.0	-
5.0	BRETE0050	9.5	-
5.5	BRETE0055	10.0	-
6.0	BRETE0060	10.5	-

# Adaptive-cuff Endotracheal Tube

*CUFF CAN EXPAND AND SHRINK AS PATIENT BREATHES*

Size: 6.0#-8.0#

- Avoid damage caused by ischemia of airway mucosa.
- Eliminate complication caused by improper cuff pressure.

Normal			
SIZE	CAT. NO	SIZE	CAT. NO
6.0	AET0060	7.5	AET0075
6.5	AET0065	8.0	AET0080
7.0	AET0070	8.5	-

Reinforced			
SIZE	CAT. NO	SIZE	CAT. NO
6.0	ARET0060	7.5	ARET0075
6.5	ARET0065	8.0	ARET0080
7.0	ARET0070	8.5	-



# Lubricious Coating Endotracheal Tube

*SUPER LUBRICATING SURFACE MODIFICATION TECHNOLOGY*

Size: Normal Type 2.0#-11.0#

Reinforced Type 2.5#-10.0#

- Use lubricious coating technology to form a lubricating film.
- Reduce the friction damage to the airway mucosa during intubation and extubation.

Normal/Reinforced					
SIZE	CAT. NO	CAT. NO (REINFORCED)	SIZE	CAT. NO	CAT. NO (REINFORCED)
2.0	LET0020	-	6.5	LET0065	LRET0065
2.5	LET0025	LRET0025	7.0	LET0070	LRET0070
3.0	LET0030	LRET0030	7.5	LET0075	LRET0075
3.5	LET0035	LRET0035	8.0	LET0080	LRET0080
4.0	LET0040	LRET0040	8.5	LET0085	LRET0085
4.5	LET0045	LRET0045	9.0	LET0090	LRET0090
5.0	LET0050	LRET0050	9.5	LET0095	LRET0095
5.5	LET0055	LRET0055	10.0	LET0100	LRET0100
6.0	LET0060	LRET0060	10.5	-	-





## TUOREN MEDICAL

- +86-373-8605444
- info@bjtuoren.com; info@etuoren.com
- Tuoren Industrial Zone, Changyuan, Henan, China
- [www.tuoren.com](http://www.tuoren.com)

