

Berman Oropharyngeal Airways

Product Data Sheet

PRODUCT: Berman Oropharyngeal Airways



DESCRIPTION / FUNCTION

Flexicare's Berman Airways are oropharyngeal airways designed to maintain the patency of a patient's airway through the oral cavity by holding forward the base of the tongue, preventing it falling backwards and obstructing the airway.

WARNINGS

- Do not use if packaging is open or damaged
- Single use device. Do not re-use, sterilise or autoclave
- After use, this product is contaminated waste and may be a potential biohazard. Handle and dispose of in accordance with hospital policy and applicable local guidelines and regulations

SPECIFICATION / KEY FEATURE

- Side channels enable the use of suction catheters without obstructing the airway, whilst allowing for additional air flow

PRODUCT RANGE & PART NUMBERS

038-90-030	Length 3cm, Size 0000
038-90-040	Length 4cm, Size 000
038-90-050	Length 5cm, Size 00
038-90-060	Length 6cm, Size 0
038-90-070	Length 7cm, Size 1
038-90-080	Length 8cm, Size 2
038-90-090	Length 9cm, Size 3
038-90-100	Length 10cm, Size 4
038-90-110	Length 11cm, Size 5
038-90-120	Length 12cm, Size 6

MATERIALS

Component	Material
Airway	High Density Polyethylene (PE)

LATEX CONTENT

Flexicare's Berman Airways do not contain natural rubber Latex.

DEHP CONTENT

Flexicare's Berman Airways do not contain phthalate DEHP.

SINGLE-USE

Flexicare's Berman Airways are single use devices.
Do not Re-use, sterilise or Autoclave.

STERILITY

Flexicare's Berman Airways are supplied non-sterile.

STORAGE

Store in a cool, dry place out of direct sunlight.

SHELF LIFE

Flexicare declares a shelf life of five years from the date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product packaging.

DISPOSAL CONSIDERATIONS

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

PACKAGING MATERIALS

Primary - Polybag - Polyethylene
Secondary – Cardboard Carton / Box, Quantity per Carton = 50