

# LarySeal™ Blue

## Laryngeal Mask Airway

Product Data Sheet

**PRODUCT:** LarySeal™ Blue



**DESCRIPTION / FUNCTION**

The LarySeal™ Blue Laryngeal Mask Airway is used for patients who require anaesthesia, artificial respiration and other types of assisted respiration. It is an alternative method of maintaining the patency of the airway against the traditional use of an Endotracheal Tube.

The LarySeal™ Blue was introduced to provide all the benefits of the reusable LarySeal Multiple whilst enhancing protection against cross contamination as it is a single use product and the first silicone single patient use laryngeal mask airway.

**WARNINGS**

- Do not use if packaging is open or damaged
- Consult Instructions for Use
- Single use. Do not reuse
- After use, this product is contaminated clinical waste and may be a potential biohazard. Handle and dispose of in accordance with hospital policy and applicable local guidelines and regulations.

## SPECIFICATION / KEY FEATURES

- Medical Grade Silicone
- Single use
- Thickness of the main tube wall reduces the risk of crushing
- The matte surface of the LarySeal™ Blue cuff provides good adhesion to the oropharyngeal area and reduces the risk of movement or dislocation of the cuff
- ISO 5356 compliant 15mm male connector provides a universal connection.
- Supplied sterile by Ethylene Oxide Gas
- Thin walled pilot balloon indicates the inflation of the cuff
- Ease of identification of sizing by colour coded pilot balloon
- Latex free
- DEHP free

## PRODUCT RANGE & PART NUMBERS

Product Code	Description	Quantity
038-94-210	LarySeal Blue Size 1 < 5kg ①	10
038-94-215	LarySeal Blue Size 1.5 5kg -10kg ①.5	10
038-94-220	LarySeal Blue Size 2 10kg - 20kg ②	10
038-94-225	LarySeal Blue Size 2.5 20kg - 30kg ②.5	10
038-94-230	LarySeal Blue Size 3 30kg - 50kg ③	10
038-94-240	LarySeal Blue Size 4 50kg - 70kg ④	10
038-94-250	LarySeal Blue Size 5 70kg - 100kg ⑤	10

### New LarySeal Blue in Slimline Pack

Product Code	Description	Quantity
038-94-610	LarySeal Blue Size 1 < 5kg ①	20
038-94-615	LarySeal Blue Size 1.5 5kg -10kg ①.5	20
038-94-620	LarySeal Blue Size 2 10kg - 20kg ②	20
038-94-625	LarySeal Blue Size 2.5 20kg - 30kg ②.5	20
038-94-630	LarySeal Blue Size 3 30kg - 50kg ③	20
038-94-640	LarySeal Blue Size 4 50kg - 70kg ④	20
038-94-650	LarySeal Blue Size 5 70kg - 100kg ⑤	20

## MATERIALS

Component	Material
Cushion	Silicone
Shell	Silicone
Main Tube	Silicone
Inflation Tube	Silicone
Pilot Balloon	Silicone
15mm Connector	Acrylonitrile Butadiene Styrene (ABS)
Valve	Polyester / Stainless Steel / Polyoxymethylene (POM) / Nitrile Rubber
Valve tag	Polypropylene (PP)

## **LATEX CONTENT**

Flexicare LarySeal™ Blue does not contain natural rubber latex.

## **DEHP CONTENT**

Flexicare LarySeal™ Blue does not contain phthalate DEHP.

## **SINGLE USE**

Flexicare LarySeal™ Blue is single use.

## **STERILITY**

Flexicare LarySeal™ Blue is supplied sterile by Ethylene Oxide (EO) gas method.

## **STORAGE**

Store in a cool, dry place out of direct sunlight.

## **SHELF LIFE**

Flexicare declares a shelf life of five years from the date of sterilisation. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product pouch.

## **DISPOSAL CONSIDERATIONS**

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

## **PACKAGING MATERIALS**

Primary - Top Web – EO Permeable High-density Polyethylene (HDPE) filaments  
Bottom Web – Polyvinyl Chloride (PVC)  
Secondary – Carton / Box Cardboard