

Flexible (McCoy) Reusable Fibre Optic Laryngoscope Blades

Product Data Sheet

PRODUCT: Flexible (McCoy) Reusable Fibre Optic Laryngoscope Blades



DESCRIPTION / FUNCTION

Flexicare's Flexible Reusable Fibre Optic Laryngoscope Blades are designed to be connected to fibre optic handles to channel the light from the source to the tip of the blade. They have a hinged tip operated by a lever mechanism. This elevates the epiglottis while reducing the amount of force required. A flexible blade has been shown to improve the view during difficult intubations and in patients wearing cervical collars.

WARNINGS

- Do not use if packaging is open or damaged.
- Consult Instructions for Use
- After use, this product is contaminated waste and may be a potential biohazard. Handle and dispose of in accordance with hospital policy and applicable local guidelines and regulations.

SPECIFICATION / KEY FEATURES

- Flexicare Fibre Optic Laryngoscope Blades are compatible with all fibre optic laryngoscope handles that comply with the standard ISO 7376
- The design matches that of traditional McCoy blades in respect of length and shape.
- Suitable for Multiple patients and can be autoclaved.
- The locking block assembly engages with all fibre optic handles.
- Light Guide provides excellent light transmission > 600 lux
- Higher outputs can be achieved depending on light source of Fibre Optic Handle.
- Rounded tip reduces the risk of soft tissue damage
- Size identification and on side of locking block
- Flexible tip controlled by lever at rear of blade

PRODUCT RANGE & PART NUMBERS

040-835	Flexible Fibre Optic Blade, Size 3
040-836	Flexible Fibre Optic Blade, Size 4

MATERIALS

Component	Material
Blade	Stainless Steel
Lever	Stainless Steel
Locking Block	Stainless Steel
Light Guide	Glass
Light Guide Sheath	Stainless Steel
Locking bearings	Stainless Steel

LATEX CONTENT

Flexicare's Flexible Fibre Optic Blades do not contain natural rubber Latex.

DEHP CONTENT

Flexicare's Flexible Fibre Optic Blades do not contain phthalate DEHP.

SINGLE-USE

Flexicare's Flexible Fibre Optic Blades are Reusable devices.

STERILITY

Flexicare's Flexible Fibre Optic Blades are supplied non-sterile.

STORAGE

Store in a cool, dry place out of direct sunlight.

SHELF LIFE

Flexicare declares a shelf life of five years from the date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product pouch.

DISPOSAL CONSIDERATIONS

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

PACKAGING MATERIALS

Primary - Foam insert within protective carry case
Secondary – Cardboard Carton / Box, Quantity per carton = 1