

Fibre Optic Reusable Laryngoscope Blades

Miller Blades
Product Data Sheet

PRODUCT: Miller Reusable Fibre Optic Laryngoscope Blades



DESCRIPTION / FUNCTION

Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades are designed to be connected to fibre optic laryngoscope handles to channel the light from the source to the tip of the blade. Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades are straight and low profile in design, with a narrower flange than other laryngoscope blades. The narrower flange is designed to reduce risk of soft tissue trauma and the curved tip facilitates easy anterior lifting of the epiglottis and greater exposure of the larynx.

Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades are supplied non-sterile.

Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades are intended for use in both pre-hospital and hospital environments.

WARNINGS

- Use care, do not point the light guide directly at the eye while in operation. The light output can cause severe eye discomfort.
- Do not use if packaging is open or damaged.
- Consult Instructions for Use
- After use, this product is contaminated waste and may be a potential biohazard. Handle in accordance with hospital policy, applicable local guidelines and regulations and manufacturer reprocessing instructions.

SPECIFICATION / KEY FEATURES

- Latex Free.
- DEHP Free.
- Fibre optic (green specification) handle.
- Hook-on fitting of handle complies with BS EN ISO 7376:2009.
- Locking block of blades complies with BS EN ISO 7376:2009.

PRODUCT RANGE & PART NUMBERS

Product Code	Description	Quantity
Miller – Reusable Fibre Optic		
040-849	Miller fibre optic reusable laryngoscope blade size 00	1
040-850	Miller fibre optic reusable laryngoscope blade size 0	1
040-851	Miller fibre optic reusable laryngoscope blade size 1	1
040-852	Miller fibre optic reusable laryngoscope blade size 2	1
040-853	Miller fibre optic reusable laryngoscope blade size 3	1
040-854	Miller fibre optic reusable laryngoscope blade size 4	1

MATERIALS

Component	Material
Blade	
Blade	Stainless steel
Locking Block	Stainless steel
Light Guide	Glass
Light Guide Sheath	Stainless steel
Locking Bearing	Stainless Steel

LATEX CONTENT

Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades do not contain natural rubber latex.

DEHP CONTENT

Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades do not contain phthalate DEHP.

MULTIPLE USE

Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades are multiple use devices, intended for multiple patient use. To be cleaned and sterilised in accordance with hospital policy, applicable local guidelines and regulations and manufacturer reprocessing instructions.

STERILITY

Flexicare's Miller Reusable Fibre Optic Laryngoscope Blades re supplied non-sterile.

STORAGE

Store in a cool, dry place out of direct sunlight.

SHELF LIFE

Flexicare declares a shelf life of five years from the date of manufacture. This is based on the stability of the devices' components and raw materials sourced.

PACKAGING MATERIALS

Primary - Carton / Box - Cardboard