

Catheter Mounts

Product Data Sheet

PRODUCT Catheter Mounts

Corrugated Tubing



Smooth Bore Tubing



Flexi-lock Tubing

DESCRIPTION / FUNCTION

Flexicare's catheter mounts are designed to provide a flexible connection between an intubation device and a breathing circuit and reduce the drag on an ET Tube or Laryngeal Mask by transferring the weight of the breathing system away from the patient. The flexibility in the tube enables the position / approach of the main breathing circuit to be selected for the convenience of clinicians and the comfort of the patient. A wide range of catheter mounts are available with Standard, Flexi-lock or Smooth Bore tubing with a variety of fixed or swivel elbows and connectors offering total flexibility.

PRE-CAUTIONS / WARNINGS

- DO NOT use if the packaging is damaged
- DO NOT reuse

SPECIFICATION / KEY FEATURES

- Conforming to BS EN 12342 Breathing tubes intended for use with Anaesthesia apparatus and Ventilators
- 15mm standard or Flexi-lock corrugated tubing allows ultimate flexibility during operative and clinical procedures when patient manoeuvring may be required
- Wide range of fixed or swivel elbows with or without ports providing total flexibility for patient airway management
- 22mm female end connectors allow for secure fitting to all Y-pieces or filters
- 15mm female end connectors allow for secure fitting to all ET tube or Laryngeal Mask Airway
- Crush resistant
- Flow is still maintained even if twisted through 180°
- Latex free
- DEHP free

PRODUCT RANGE & PART NUMBERS

Standard 15mm Corrugated Tubing

038-61-301	Standard 15mm Tube Catheter Mount, 22mm F & 15mm F Connectors
038-61-302	Standard 15mm Tube Catheter Mount with Low Profile Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-303	Standard 15mm Tube Catheter Mount with Luer Lock Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-304	Standard 15mm Tube Catheter Mount with Low Profile Elbow, Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors
038-61-305	Standard 15mm Tube Catheter Mount with Double Swivel Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-306	Standard 15mm Tube Catheter Mount with Double Swivel Elbow, Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors
038-61-301/640	Standard 15mm Tube Catheter Mount with Detachable Low Profile Elbow, 22mm F & 15mm F/22mm Connectors
038-61-301/641	Standard 15mm Tube Catheter Mount with Detachable Luer Lock Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-301/642	Standard 15mm Tube Catheter Mount with Detachable Low Profile Elbow, Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors
038-61-301/643	Standard 15mm Tube Catheter Mount with Low Profile Elbow, Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors
038-61-301/645	Standard 15mm Tube Catheter Mount with Double Swivel Elbow, 22mm F & 15mm F/22mm M Connectors

Smooth Bore Tubing

038-61-321	Smooth Bore Catheter Mount, 22mm F & 15mm F Connectors
038-61-322	Smooth Bore Catheter Mount with Low Profile Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-323	Smooth Bore Catheter Mount with Luer Lock Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-324	Smooth Bore Catheter Mount with Elbow, Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors
038-61-325	Smooth Bore Catheter Mount with Double Swivel Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-326	Smooth Bore Catheter Mount with Double Swivel Elbow, Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors

Flexi-lock Tubing

038-61-310	Flexi-lock 15mm Tube Catheter Mount, 15mm M & 15mm F Connector
038-61-311	Flexi-lock 22mm Tube Catheter Mount 22mm F & 15mm F Connector
038-61-312	Flexi-lock Catheter Mount with Low Profile Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-313	Flexi-lock Catheter Mount with Luer Lock Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-314	Flexi-lock Catheter Mount with Elbow and Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors
038-61-315	Flexi-lock Catheter Mount with Double Swivel Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-316	Flexi-lock Catheter Mount with Double Swivel Elbow Elastomeric Cap, 22mm F & 15mm F/22mm Connectors
038-61-310/640	Flexi-lock Catheter Mount with Detachable Low Profile Elbow, 22mm F & 15mm F/22mm M Connectors
038-61-310/641	Flexi-lock Catheter Mount with Detachable Luer Lock Elbow, 22mm F & 15mm M Connector
038-61-310/642	Flexi-lock Catheter Mount with Detachable Elbow, Elastomeric Cap, 22mm F & 15mm F/15mm M Connector
038-61-310/643	Flexi-lock Catheter Mount with Detachable Double Swivel Elbow, 22mm F & 15mm F/22mm M Connector

038-61-310/645	Flexi-lock Catheter Mount with Detachable Double Swivel Elbow, Elastomeric Cap, 22mm F & 15mm F/22mm M Connectors
----------------	---

MATERIALS

Description	Material
Corrugated Tubing	Polyethylene (PE)
Smooth Bore Tubing	Polyester Elastomer
Flexi-lock Tubing	Thermoplastic Elastomer
15m/22F Connectors	Acrylonitrile Butadiene Styrene (ABS)
Elastomeric Cap	Thermoplastic Elastomer
Low Profile Elbow	Polypropylene (PP)
Luer Lock Cap & Strap	Polypropylene (PP)
Double Swivel Elbow Body	Polyethylene (PE)
Double Swivel Connector	Acrylonitrile Butadiene Styrene (ABS)

LATEX CONTENT

Flexicare Catheter Mounts do not contain natural rubber Latex.

DEHP CONTENT

Flexicare Catheter Mounts do not contain phthalate DEHP.

SINGLE-USE

Flexicare Catheter Mounts are single use.

STERILISATION

Flexicare Catheter Mounts are supplied non-sterile.

STORAGE CONDITIONS

Store in a cool, dry place out of direct sunlight.

SHELF LIFE

Flexicare declares a shelf life of five years from date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product packaging.

DISPOSAL CONSIDERATIONS

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

PACKAGING

Primary – Polybag Polyethylene (PE)
 Secondary – Carton Cardboard, Quantity per carton = 50