


Suction Tubing

Product Data Sheet

PRODUCT Suction Tubing**DESCRIPTION / FUNCTION**

Flexicare's Suction Tubing is designed to provide a flexible conduit to connect receiver and vacuum source to device at patient end used to remove unwanted fluids during medical procedures and treatments. The suction tube is resistant to collapse under vacuum.

PRE-CAUTIONS / WARNINGS

- DO NOT use if the packaging is damaged
- DO NOT re-sterilise

SPECIFICATION / KEY FEATURES

- Available with internal diameter (ID) from 5mm to 7mm and lengths from 2m through to 30m
- Wall thickness prevents the tube from collapsing when under high pressure
- Easy to open sterile outer peel pouch
- Universal female connector colour coded for ID size
- Spare male to male connector for connection to female connectors, or to connect two tubes if extended length is required
- Latex free
- DEHP free

PRODUCT RANGE & PART NUMBERS

1020-25	7mm ID Suction Tube Sterile, 2m with Yellow F/F Connector	●
1020-6	6mm ID Suction Tube Sterile, 2m with Green F/F Connector	●
1024-25	7mm ID Suction Tube Sterile, 2.4m with Yellow F/F Connector	●
1027-25	7mm ID Suction Tube Sterile, 2.7m with Yellow F/F Connector	●
1030-25	7mm ID Suction Tube Sterile, 3m with Yellow F/F Connector	●
1030-5	5mm ID Suction Tube Sterile, 3m with Blue F/F Connector	●
1030-6	6mm ID Suction Tube Sterile, 3m with Green F/F Connector	●
1035-25	7mm ID Suction Tube Sterile, 3.5m with Yellow F/F Connector	●
1037-25	7mm ID Suction Tube Sterile, 3.7m with Yellow F/F Connector	●
10830-5	Suction Tube Non-sterile 30m Coil, 5mm ID with Bubble Non Conductive	
10830-7	Suction Tube Non-sterile 30m Coil, 5mm ID with Bubble Non Conductive	

MATERIALS

Tubing	Polyvinyl Chloride (PVC)
Connectors	Polyvinyl Chloride (PVC)

LATEX CONTENT

Flexicare's Suction Tubing does not contain natural rubber Latex.

DEHP CONTENT

Flexicare's Suction Tubing does not contain phthalate DEHP.

SINGLE-USE

Flexicare's Suction Tubing is single use.

STERILISATION

Flexicare's Suction Tubing is supplied sterile by Ethylene Oxide gas.
30m Suction Tube coils are supplied non-sterile.

STORAGE CONDITIONS

Store in a cool, dry place out of direct sunlight.

SHELF LIFE

Flexicare declares a shelf life of five years from date of sterilisation. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product packaging.

DISPOSAL CONSIDERATIONS

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

PACKAGING

Primary – Sterile Pouch or Non-sterile Polybag
Secondary – Carton Cardboard

BOX / CARTON

Product	Dimensions (mm)	Qty per carton
1020-25	W 385mm L 190mm H 315mm	25
1020-6	W 385mm L 190mm H 315mm	25
1024-25	W 385mm L 190mm H 315mm	25
1027-25	W 385mm L 190mm H 315mm	25
1030-25	W 385mm L 190mm H 315mm	25
1030-5	W 385mm L 190mm H 315mm	25
1030-6	W 385mm L 190mm H 315mm	25
1035-25	W 385mm L 190mm H 315mm	25
1037-25	W 385mm L 190mm H 315mm	25
10830-5	W 312mm L 130mm H 305mm	1
10830-7	W 312mm L 130mm H 305mm	1