


Superior Safety Design

- Low risk of accidental spillage of liquid or contact
- DEHP free & Biologically decomposable




Hydrophobic Filter

- Serve as a combined bacterial filter and overflow protection
- High Bacterial Aerosol Filtration Efficiency

Feature	BFE*	Air Flow
Self-sealing	>99.98%	29.2 scfh

* Bacterial Filtration Efficiency
(Test method conform to ASTM F2102)



Overflow protection

- Automatically shut off when waste fluid get full in liner
- Avoid waste from being inhaled to medical suction system

Avoid premature shut-off

- High air flow rate