



OPERATING ROOM GEL POSITIONING PADS

ORIGINAL BLUE GEL POSITIONER BROCHURE & CATALOGUE



HEAD / NECK / FACE_



OPERATING ROOM GEL POSITIONING PADS

Closed Head Ring.

Protects the head, face, neck and ears during all type of surgery.



- 1) **AP001** / 9x4x2 cm / 3.5x1.6x0.8 inch
- 2) **AP002** / 14x6x3.5 cm / 5.5x2.4x1.4 inch
- 3) **AP003** / 20x7.5x4.5 cm / 7.9x3.0x1.8 inch
- 4) **AP004** / 20x7.5x7.5 cm / 7.9x3.0x3.0 inch

Open Head Ring.

Designed to support a patient's head and gace during surgery with easy access for surgical intubation.



- 1) **AP011** / 9x4x2 cm / 3.5x1.6x0.8 inch
- 2) **AP012** / 14x6x3.5 cm / 5.5x2.4x1.4 inch
- 3) **AP013** / 20x7.5x4.5 cm / 7.9x3.0x1.8 inch
- 4) **AP014** / 20x7.5x7.5 cm / 7.9x3.0x3.0 inch

Open Head Ring with Cut.

The padless open head ring provides flexibility ti accommodate various patients. The soft gel material adjusts to any patient's head and offers sufficient support and comfort during operations.




- 1) **AP011/C** / 9x4x2 cm / 3.5x1.6x0.8 inch
- 2) **AP012/C** / 14x6x3.5 cm / 5.5x2.4x1.4 inch
- 3) **AP013/C** / 20x7.5x4.5 cm / 7.9x3.0x1.8 inch



Prone Head Rest.

Unique design provides excellent comfort, stability and proper support during surgery.



- 1) **AP022** / 28x23x14.5 cm / 11x9.1x5.7 inch
- 2) **AP023** / 23x19x13 cm / 9.1x7.5x5.1 inch



Head Pad.

Head pad protects the head, face, neck and ears also provides comforts and stable positioning during all kinds of surgery.




1) AP110 / 27x23x7 cm / 10.6x9.1x2.8 inch



Contoured Head Rest.

This pad is designed for procedures in the supine position requiring optimum stability of the patient head. The soft gel material provides even weight distribution to protect the patient from developing pressure sores.




1) AP081 / 19x19x6 cm / 7.5x7.5x2.4 inch



Supine Head Rest.

Designed for patient in supine position to relieve pressure on ears, neck and back of the head.



1) AP112 / 27x24x9 cm / 10.6x9.4x3.5 inch



Face Pad.

This versatile pad can be used in various procedures where endotracheal tubes are employed. This pad offers maximum protection to the patient forehead from pressure during long hours of operation.





1) AP021 / 16x6x1.6 cm / 6.3x2.4x0.6inch

ARM / ELBOW_



Armboard Pad.

Armboard pad provides excellent pressure relief during various surgical procedures.



1) **AP071** / 40x12.5x2 cm / 16x4.9x0.8 inch

2) **AP072** / 60x12.5x2 cm / 24x4.9x0.8 inch



Contoured Armboard Pad with Foam.

Clearview's contoured armboard pads are designed to provide added support and stability.



1) **AP073** / 40x12.5x6 cm / 16x5x2.4 inch

2) **AP074** / 60x12.5x6 cm / 24x5x2.4 inch


3) **AP075** / 50x12.5x6 cm / 20x5x2.4 inch

4) **AP076** / 63x12.5x6 cm / 24.8x5x2.4 inch



Arm Protector.

Arm protector provides stable pressure relief with adjustable velcro.



1) **AP310** / 38x12x0.8 cm / 15x4.7x0.3 inch

2) **AP311** / 46x21x0.8 cm / 15.1x8.3x0.3 inch



Gel Fastening Strap.


The gel fastening strap provides effective protection for elbows and ankles through the use of soft and comfortable polyurethane gel.



1) **FB01** / 33x7x1.5 cm / 13x2.75x0.6 inch

2) **FB02** / 36x7x1.5 cm / 14.2x2.75x0.6 inch

Hand & Wrist Supporter



1) **AP005** / 28x8.3x6 cm / 11x3.3x2.4 inch

Contoured Armboard Pad.



1) **AP077** / 52x15x5.5 cm / 20.5x6x2.1 inch



Heel Support Pad.

Light weight heel pads provide maximum comfort and support for heels.
*2pcs/1set



1) AP031 / 20x11x7 cm / 8x4.3x2.75 inch



Heel Cups without Strap.

Delivers maximum heel protection for presurgery and postsurgery applications.
*2pcs/1set



1) AP032 / 17x11x5 cm / 6.7x4.1x2 inch



Comfort Heel Cups with Strap.

This product provides support for the heels especially during recovery.
*2pcs/1set



1) AP033 / 17x18x1.3 cm / 6.7x7x0.5 inch



Comfort Heel Cups with Strap.

This product provides support for the heels especially during recovery.
*2pcs/1set



1) AP034 / 17x23x1.3 cm / 6.7x9x0.5 inch

ANKLE / HEEL / FOOT_



Ankle Strap Stirrup.

Ankle strap stirrups provide protection for a patient's skin for procedures in the lithotomic position.



1) AP904 / 31x4x1 cm / 12.2x1.6x0.4 inch

2) AP905 / 44x4.5x1 cm / 17.5x1.8x0.4 inch

BODY_



Chest Rolls(Adult).

Virtually the same as pediatric chest rolls with the same function to provide pressure relief and pressure sore prevention but in bigger sizes to accommodate adult patients.



1) AP041 / 30x15x10 cm / 11.8x5.9x3.9 inch

2) AP042 / 40x15x10 cm / 15.7x5.9x3.9 inch

3) AP043 / 50x15x10 cm / 19.7x5.9x3.9 inch

4) AP044 / 30x10x8.5 cm / 11.8x3.9x3.3 inch

5) AP045 / 40x10x8.5 cm / 15.7x3.9x3.3 inch

6) AP046 / 50x10x8.5 cm / 19.7x3.9x3.3 inch

7) AP047 / 50x18x15 cm / 19.7x7.1x5.9 inch



Chest Rolls(Pediatric).

Virtually the same as pediatric chest rolls with the same function to provide pressure relief and pressure sore prevention.



1) AP047-1 / 15x5x5 cm / 5.9x2x2 inch

2) AP047-2 / 20x5x5 cm / 7.9x2x2 inch

3) AP047-3 / 25x5x5 cm / 9.8x2x2 inch

4) AP047-4 / 30x5x5 cm / 11.8x2x2 inch

Chest Rolls(Pediatric).



1) AP048-1 / 15x7.5x7.5 cm / 5.9x3x3 inch

2) AP048-2 / 20x7.5x7.5 cm / 7.9x3x3 inch

3) AP048-3 / 25x7.5x7.5 cm / 9.8x3x3 inch

4) AP048-4 / 30x7.5x7.5 cm / 11.8x3x3 inch



Contoured Chest Rolls.

The chest roll has with flat bottom provides stability and great utility.



- 1) AP049-1 / 30x16x13/9 cm / 11.8x6.3x5.1/3.5 inch
- 2) AP049-2 / 40x16x13/9 cm / 15.7x6.3x5.1/3.5 inch
- 3) AP049-3 / 50x16x13/9 cm / 19.7x6.3x5.1/3.5 inch



Universal Patient Positioner.

Universal chest positioners provide positioning aid during various surgical procedures. These pressure-distributive pads protect patients against bedsores.



- 1) AP091 / 16x8x4.5 cm / 6.3x3.1x1.7 inch
- 2) AP092 / 33x7x5 cm / 13x2.7x2.0 inch
- 3) AP093 / 33x13x6 cm / 13x5x2.4 inch
- 4) AP094 / 36x12.5x7 cm / 14.2x4.9x2.75 inch
- 5) AP095 / 25x12.5x5 cm / 9.8x4.9x2 inch
- 6) AP096 / 38x12.5x6.3 cm / 15x4.9x2.5 inch
- 7) AP097 / 50x12.5x7 cm / 20x4.9x2.75 inch



Dome Positioner.

The dome positioner is designed for use during surgery to protect the patient while in the prone or lateral position.




- 1) AP051 / 35x15x7.53 cm / 13.8x6x3 inch



Graduated Thoracic Positioner.

This positioner with its sloping design is ideal for procedures requiring full chest hyperextension. This pad provide pressure relief and nerve protection.




- 1) AP052 / 23x36x3.8 cm / 9x14.2x1.5 inch

BODY_



Fracture Table Pad.

This pad, when being used with Fracture Table post pad, effectively reduces pressure and shear at the sacral area.



1) AP054 / 14x20x1.3 cm / 5.5x8x0.5 inch



Lateral Positioner.

It is designed to accommodate both the shoulder and upper arm and provides exceptional pressure relief and nerve protection.



1) AP101 / 71x50x10 cm / 31x20x3.9 inch



Chest Positioning Set.

Chest Positioner have a contoured design patterned to "uniformly conform to a patient's shoulder-to-chest area" with round edges for the axillar. They provide optimum upper body positioning while alleviating pressure on the thorax area.



1) AP050-1 / 71/58x20x7.6 cm / 28x7.9x3 inch

2) AP050-2 / 71/58x20x7.6 cm / 28x7.9x3 inch



Patient Positioning Belt.


Positioning belt provides safety, stability and tissue protection during surgery.



1) AP601 / 51x9x0.8 cm / 20x3.3x0.3 inch

Slanting Shaped Pad.

The pad is designed to support and protect lumbar area, hepatic area and splenic area in lateral supine surgeries or supine surgeries.



1) AP102-3 / 30x22x5/1 cm / 12x9x2/0.3 inch


2) AP102-6 / 50x12x10/4 cm / 20x5x4/1.5 inch

Sacral Protector.

The sacral area is an area of serious risk if left unprotected. The protector has been designed to provide support in the supine and lithotomy positions by relieving pressure and focusing on peripheral support.



1) AP103-1 / 40x38x2.5 cm / 15.7x15x1 inch




Relton Frame Pad Set.

This 4-piece set is designed to maintain the patient spinal position during surgery. These pads provides exceptional protection against pressure sores and nerve or tissue damage for patinents undergoing long surgical procedures. [*4pcs/1se](#)



1) AP111 / 19x17.5x10.5 cm / 7.5x6.9x4.13 inch



Wilson Frame Pad.

Wilson Frame Pads help to maintain a patient's flexible spine position.



1) AP210 / 75x26x2.5 cm / 30x10x1 inch

2) AP211 / 75x26x1 cm / 30x10x0.4 inch



Wilson Frame Set.

The Wilson Frame is designed to provide the same support that the matal frames dowwithout the bulky weight and high cost. The Wilson Gel Pad fits on the Wilson Grame Pad perfectly, and is secured on the frame with velcro and straps.



1) AP212 / 72.5x17.5x2 cm / 28.5x6.9x0.8 inch - 2 ea
(Gel Pad) 18.5x12x2 cm / 7.3x4.7x0.8 inch
18.5x12x2 cm / 7.3x4.7x0.8 inch - w/c

2) AP213 / 60x49x15.5 cm / 23.6x19.3x6.1 inch
(Frame)



Prone Head Rest.


Prone Head Rest is designed so that the patient's head can comfortably rest face down with respiratory tubes in place. The unique design provides excellent support, stability, comfort and protection during surgery in the prone position. with opening both sides.



1) AP020 / 15x28x15 cm / 5.9x11x5.9 inch

WEIGHT : 5 kg
MATERIAL : Aluminum, PU, Gel Pad
FEATURES : Width Adjustable


WILSON FRAME_



Wilson Frame.

Wilson Frame Pads are designed to provide pressure relief for patients positioned on a Wilson Frame, used to maintain patients in a flexed spine position.

*Height and Width are adjustable as user's body.



1) AP214 / 40x69xx21~26.5 cm / 16x27x8.2~10.4 inch


WEIGHT : 19.4 kg

MATERIAL : Aluminum, PU, Gel Pad


FEATURES : Available X-Ray measure
Height and Width adjusted

HOW TO USE.

1. Please insert the controller to crew as below image.



2. Turn the controller.




*Width adjustment

:To make it close,
please turn the controller to the right.

*Height adjustment

:To make it lower,
please turn the controller to the right.





*Width adjustment

:To widen the width,
please turn the controller to the left.

*Height adjustment


:To make it higher,
please turn the controller to the left.





Fracture Table Hip Pad.

It is designed for use on the fracture table plate to reduce pressures and shear at the sacral area.



- 1) AP055 / 18x18x1.3 cm / 7x7x0.5 inch
- 2) AP056 / 33x30x0.8 cm / 13x11x0.3 inch



Fracture Table Post Pad.

Designed to be used with Fracture Table Hip Pad, this pad helps reduce pressure and shear at the sacral area.



- 1) AP556 / 15x13x2.5 cm / 6x5x1 inch



Crutch Stirrup Pad.

Crutch stirrup pads are designed to fit both Lloyd devices and various types of stirrups.




- 1) AP801 / 31x23x1 cm / 12x9x0.4 inch
- 2) AP802 / 40x25x1 cm / 16x9.8x0.4 inch
- 3) AP803 / 52x30x1 cm / 20.4x11.8x0.4 inch
- 4) AP804 / 55x52x1 cm / 21.6x20.4x0.4 inch



Operating Table Pad with Foam.

It is possible to use in various position with light weight polyurethane foam base.



- 1) AP419 / 70x20x8 cm / 27.6x7.9x3.1 inch
- 2) AP420 / 30x15x10 cm / 11.8x5.9x3.9 inch
- 3) AP421 / 50x15x10 cm / 19.7x5.9x3.9 inch
- 4) AP422 / 60x15x10 cm / 23.6x5.9x3.9 inch

SPECIALTY PAD_



Boot Stirrup Pad.

Boot stirrup pads provide excellent protection over calf and heel area during surgery in the lithotomin position.



1) AP901 / 71x41T/30x1 cm / 28x16/12x1/2 inch



Traction Boot Pad Cover.

Used in conjunction with gel pads to protect a patient's heel.

1) AP902 / 37x24x0.5 cm / 14.5x9.5x0.2 inch




Lithotomy Pole Cover.

Placed around lithotomy poles, these pads help protect a patient's skin.




1) AP903 / 71x12x0.8 cm / 28x4.7x0.3 inch



Ez Feel Cushion

The Ez Feel Cushion is made with a high-quality polyurethane gel. The surface is solid, smooth and helps dissipate heat to provide excellent pressure distribution and comfort. The true gel material will not migrate, leak, evaporate, harden or dry out if the cushion is accidentally punctured.


- 1) **NC1616** / 40x40x1.6 cm / 16x16x0.6 inch
- 2) **NC1816** / 46x40x1.6 cm / 18x16x0.6 inch
- 3) **NC1818** / 46x46x1.6 cm / 18x18x0.6 inch
- 4) **NC2018** / 50x46x1.6 cm / 20x18x0.6 inch



Comfy Cushion

The Comfy Cushion is made from an advanced polymer gel and features an embossed design to help dissipate heat while providing outstanding pressure relief, comfort and long lasting durability. The true gel material will not migrate, leak, evaporate, harden, or dry out if the cushion is accidentally punctured. The material is highly resistant to bacteria and bodily fluids and can be easily cleaned with warm soapy water or disinfectant.

- 1) **NP1616** / 40x40x1 cm / 16x16x0.4 inch
- 2) **NP1816** / 46x40x1 cm / 18x16x0.4 inch
- 3) **NP1818** / 46x46x1 cm / 18x18x0.4 inch
- 4) **NP2018** / 50x46x1 cm / 20x18x0.4 inch




Ultra Cushion

The new Ultra / Air design is a revolutionary new technology for wheelchair cushions. The quilted design helps dissipate heat while providing outstanding pressure relief, comfort and long lasting durability. The true gel material will not migrate, leak, evaporate, harden, or dry out if the cushion is accidentally punctured. The material is highly resistant to bacteria and bodily fluids and can be easily cleaned with warm soapy water or disinfectant.

- 1) **NAL1616** / 40x40x1 cm / 16x16x0.4 inch
- 2) **NAL1816** / 46x40x1 cm / 18x16x0.4 inch
- 3) **NAL1818** / 46x46x1 cm / 18x18x0.4 inch
- 4) **NAL2018** / 50x46x1 cm / 20x18x0.4 inch

1) Operating Table Pad



2) VAC-PAC Cover



3) Operating Table Pad Set(3 Pcs)



4) Embossed Operating Table Pad



TABLE PAD_

Table pads provide protection over the entire surface area of the operating table and are recommended for use with patients in all positions.

The pads are flexible, soft, comfortable, versatile and distribute weight evenly over the entire surface area.

Operating Table Pad.

AP301 / 23x50x1 cm / 8.5x20x0.4 inch
AP302 / 99x50x1 cm / 39x20x0.4 inch
AP303 / 50x50x1 cm / 20x20x0.4 inch
AP303-1.5 / 50x50x1.5 cm / 20x20x0.6 inch

AP304 / 50x25x1 cm / 20x10x0.4 inch
AP304-1.5 / 50x25x1.5 cm / 20x10x0.6 inch
AP305 / 10x38x1 cm / 4x15x0.4 inch
AP306 / 42x15x1 cm / 16.5x6x0.4 inch

AP307 / 10x61x1 cm / 4x24x0.4 inch
AP308 / 45x50x1 cm / 18x20x0.4 inch
AP308-1.5 / 45x50x1.5 cm / 18x20x0.6 inch
AP402 / 183x50x1 cm / 72x20x0.4 inch

AP402-1.5 / 183x50x1.5 cm / 72x20x0.6 inch
AP401 / 193x86x1 cm / 76x33.8x0.4 inch
AP403 / 117x50x1 cm / 46x20x0.4 inch
AP403-1.5 / 117x50x1.5 cm / 46x20x0.6 inch

AP404 / 160x43x1 cm / 63x17x0.4 inch
AP405-1.5 / 15x35x1.5 cm / 6x14x0.6 inch
AP406-1.5 / 70x50x1.5 cm / 27.5x20x0.6 inch
AP407-1.5 / 70x45x1.5 cm / 27.5x18x0.6 inch

AP408-1.5 / 15x50x1.5 cm / 6x20x0.6 inch
AP053 / 31x66x1.3 cm / 12x26x0.5 inch
AP555 / 15x66x1.3 cm / 6x26x0.5 inch
AP409W/C / 117x50x1 cm / 46x20x0.4 inch

AP409W/C-1.5 / 117x50x1.5 cm / 46x20x0.6 inch
AP410W/C / 100x50x1 cm / 40x20x0.4 inch
AP410W/C-1.5 / 100x50x1.5 cm / 40x20x0.6 inch
AP411W/C / 50x50x1 cm / 20x20x0.4 inch

AP411W/C-1.5 / 50x50x1.5 cm / 20x20x0.6 inch
AP412 / 33x46x0.5 cm / 13x18x0.2 inch
AP412-1 / 33x46x1 cm / 13x18x0.4 inch
AP413 / 22x17x1 cm / 8.7x6.7x0.4 inch

AP414 / 60x60x2.5 cm / 24x24x1 inch
AP415 / 90x50x2.5 cm / 36x20x1 inch
AP416 / 75x25x2 cm / 30x10x0.8 inch
AP417 / 10x15x2 cm / 4x6x0.8 inch

AP418 / 16x8x0.6 cm / 10.7x3.2x0.2 inch
AP423 / 70x50x1.6 cm / 27.5x20x0.6 inch
AP424 / 100x50x1.6 cm / 40x20x0.6 inch

Operating Table Pad Set(3pcs).

(HEAD)

AP304 / 50x25x1 cm / 20x10x0.4 inch
AP304-1.5 / 50x25x1.5 cm / 20x10x0.6 inch

(BODY)

AP410W/C / 100x50x1 cm / 40x20x0.4 inch
AP410W/C-1.5 / 100x50x1.5 cm / 40x20x0.6 inch

(FOOT)

AP411W/C / 50x50x1 cm / 20x20x0.4 inch
AP411W/C-1.5 / 50x50x1.5 cm / 20x20x0.6 inch

Embossed Operating Table Pad.

AP061 / 50x50x1.6 cm / 20x20x0.63 inch
AP062 / 115x50x1.6 cm / 45x20x0.63 inch

AP063 / 40x23x1.6 cm / 16x9x0.63 inch
AP064 / 91x45x1.6 cm / 36x18x0.63 inch

VAC-PAC Cover.

AP701 / 97x64x1 cm / 38x25.2x0.4 inch

[illegible]



Libertech

1352, Poseungjangan-ro, Jangan-myeon,
Hwaseong-si, Gyeonggi-do, 18584 Korea

Tel : 82-70-4680-0707

Fax : 82-70-8677-8550

17-1030

