

Specification Medijet[®] 1000



Technical Specification:

- Material: Terlax, Elastollan, PP, TPE-S
- Connector Flow input:
 - 10mm OD (female) / 7,4mm ID (male),
 - material: Makrolon
- Connector pressure line:
 - Female Luer connector
 - Material PC
- Production technique: injection moulding

Connections:

- Flow input
- Pressure output
- Drug nebulisation input (6,0 mm ID)
- Flow output for Prong / mask
- Expiration opening

Dimensions:

- Total length: 17cm
- Weight: 9g

nCPAP Parameter:

- nCPAP Pressure: $\geq 4,5\text{cmH}_2\text{O}$ at 7 l/min flow