

## Specification Miniflow<sup>®</sup> 4000



### Technical Specification:

- Material: Terlux, TPE-S
- Connector Flow input/ outputs
  - 10mm OD (female) / 7,4mm ID (male)
  - material: Makrolon
- Production technique: injection moulding

### Connections:

- Flow input
- Pressure output
- Flow output for Prong / mask
- Expiration

### Dimensions:

- Total length: 16cm
- Weight: 10g

### nCPAP Parameter:

- CPAP Pressure: 0 - 20cmH<sub>2</sub>O