



# PM-2000

## Central Monitoring System

# PM-2000 Central Monitoring System

## Alarm Review



## Patient Info. Review



## Trend Data Review



## Waveform Review



### Features:

- Maximum connects 64 bedside units with one central monitoring system;
- Store & review 240-hour trend data for each bedside monitor;
- Store & review 3000 alarm events for each bedside monitor;
- Up to 20000 patient information records storage;
- Store & review 72-hour data for 64 channels full disclosure waveform;
- Various display interfaces including full screen multi-bed, half screen multi-bed and half screen key beds etc;
- Support 32 kinds of screen layout configuration;
- Waveform display: 2-CH ECG, 1-CH RESP, 1-CH PR, 1-CH EtCO2, 4-CH BISxp;
- Parameter display: ECG, RESP, ST segment, PVCs, SpO2, PR, NIBP, TEMP, EtCO2 and CSI;
- Real-time bidirectional control between bedside units and central monitoring system;
- Support single screen and dual-screen display;
- Audible and visual alarm with character alarm on the screen.

### Hardware requirement:

CPU	-----	Intel dual cores / 3GHz or above
Memory	-----	2G or above
Hard disk	-----	1T or above
Network interface card	-----	100Mb or above
Display	-----	21 inches screen or above
Graphics card	-----	NVIDIA 1G discrete graphics or above

