

# LarySeal™ Multiple

## Laryngeal Mask Airway

Product Data Sheet

## PRODUCT: LarySeal™ Multiple



### DESCRIPTION / FUNCTION

The LarySeal™ Multiple Laryngeal Mask Airway is used for patients who require anaesthesia, artificial respiration and other types of assisted respiration. It is an alternative method of maintaining the patency of the airway against the traditional use of an Endotracheal Tube.

Flexicare's silicone reusable laryngeal mask airways are designed to be sterilised up to a recommended 40 times. The LarySeal™ Multiple cuff adapts to the contour of the oropharyngeal area to provide a secure seal.

### WARNINGS

- Do not use if packaging is open or damaged
- Consult Instructions for Use
- Handle and dispose of in accordance with hospital policy and applicable local guidelines and regulations

## SPECIFICATION / KEY FEATURES

- Medical Grade Silicone
- Thickness of main tube wall reduces the risk of crushing
- Silicone allows the cuff of the mask to form itself to the contour of the oropharyngeal area providing a secure seal for airway management
- ISO 5356 compliant 15mm male connector provides a universal connection
- Durable high temperature resistant material. Recommended for up to 40 times autoclaving
- Thin walled pilot balloon indicates the inflation of the cuff
- Ease of identification of sizing by colour coded pilot balloon
- Latex free
- DEHP free

## PRODUCT RANGE & PART NUMBERS

Product Code	Description		Quantity
038-94-110	LarySeal Multiple Size 1	< 5kg	10
038-94-115	LarySeal Multiple Size 1.5	5kg -10kg	10
038-94-120	LarySeal Multiple Size 2	10kg - 20kg	10
038-94-125	LarySeal Multiple Size 2.5	20kg - 30kg	10
038-94-130	LarySeal Multiple Size 3	30kg - 50kg	10
038-94-140	LarySeal Multiple Size 4	50kg - 70kg	10
038-94-150	LarySeal Multiple Size 5	70kg - 100kg	10

## MATERIALS

Component	Material
Cushion	Silicone
Shell	Silicone
Main Tube	Silicone
Inflation Tube	Silicone
Pilot Balloon	Silicone
15mm Connector	Polycarbonate (PC)
Valve	Polypropylene (PP) / Polyoxymethylene (POM) / Nitrile Rubber / Stainless Steel
Valve Tag	Polypropylene (PP)

## **LATEX CONTENT**

Flexicare LarySeal™ Multiple does not contain natural rubber Latex.

## **DEHP CONTENT**

Flexicare LarySeal™ Multiple does not contain phthalate DEHP.

## **REUSABLE**

Flexicare LarySeal™ Multiple is reusable up to a recommended 40 times.

## **STERILITY**

Flexicare LarySeal™ Multiple is supplied sterile by Ethylene Oxide (EO) gas method.

## **STORAGE**

Store in a cool, dry place out of direct sunlight.

## **SHELF LIFE**

Flexicare declares a shelf life of five years from the date of sterilisation. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product pouch.

## **DISPOSAL CONSIDERATIONS**

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

## **PACKAGING MATERIALS**

Primary - Top Web – EO permeable High-density Polyethylene (HDPE) filaments  
Thermoform Tray – Polyvinyl Chloride (PVC)  
Secondary - Carton / Box Cardboard