

Star Lumen Oxygen Tubing

Product Data Sheet

PRODUCT Star Lumen Oxygen Tubing**DESCRIPTION / FUNCTION**

Flexicare's Star Lumen Tubing is designed to provide a low pressure tube connection between an oxygen or compressed air source and a medical device in order to provide oxygen for inhalation via the mouth and nose.

PRE-CAUTIONS / WARNINGS

- DO NOT use if the packaging is damaged
- DO NOT reuse

SPECIFICATION / KEY FEATURES

- Conforming to BS EN 13544-2 Respiratory therapy equipment. Part 2 Tubing and connectors
- Anti kink design to prevent accidental blockage of the oxygen flow
- Latex free
- DEHP free

PRODUCT RANGE & PART NUMBERS

032-10-060	Star Lumen Oxygen Tubing 2.1m with Connectors
032-10-063	Star Lumen Oxygen Tubing 2.1m with one narrow and one wide connector
032-10-059	Star Lumen Oxygen Tubing 100m coil

MATERIALS

Tubing	Polyvinyl Chloride (PVC)
Connectors	Polyvinyl Chloride (PVC)

LATEX CONTENT

Flexicare's Star Lumen Tubing does not contain natural rubber Latex.

DEHP CONTENT

Flexicare's Star Lumen Tubing does not contain phthalate DEHP.

SINGLE-USE

Flexicare's Star Lumen Tubing is single use.

STERILISATION

Flexicare's Star Lumen Tubing is supplied non-sterile.

STORAGE CONDITIONS

Store in a cool, dry place out of direct sunlight.

SHELF LIFE

Flexicare declares a shelf life of five years from date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product packaging.

DISPOSAL CONSIDERATIONS

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

PACKAGING

Primary – Polybag Polyethylene (PE)

Secondary – Carton Cardboard

BOX / CARTON

Product	Dimensions (mm)	Qty per carton
032-10-060	W 350mm L 300mm H 350mm	100
032-10-063	W 350mm L 300mm H 350mm	100
032-10-059	W 336mm L 310mm H 92mm	1