


# Oxygen Bubble Tubing

Product Data Sheet

**PRODUCT** Oxygen Bubble Tubing**DESCRIPTION / FUNCTION**

Flexicare's Oxygen Tubing is designed to provide a low pressure tube connection between an oxygen or compressed air source and a medical device in order to provide oxygen for inhalation via the mouth and nose.

**PRE-CAUTIONS / WARNINGS**

- DO NOT use if the packaging is damaged
- DO NOT reuse

**SPECIFICATION / KEY FEATURES**

- Conforming to BS EN 13544-2 Respiratory therapy equipment. Part 2 Tubing and connectors
- Anti kink design to prevent accidental blockage of the oxygen flow
- Latex free
- DEHP free

**PRODUCT RANGE & PART NUMBERS**

032-10-069	Oxygen Bubble Tubing 1.8m
032-10-067	Oxygen Bubble Tubing (3mm) 30m coil in dispensing box
032-10-068	Oxygen Bubble Tubing (3mm) 100m coil in dispensing box

**MATERIALS**

Tubing	Polyvinyl Chloride (PVC)
--------	--------------------------

### LATEX CONTENT

Flexicare's Oxygen Bubble Tubing does not contain natural rubber Latex.

### DEHP CONTENT

Flexicare's Oxygen Bubble Tubing does not contain phthalate DEHP.

### SINGLE-USE

Flexicare's Oxygen Bubble Tubing is single use.

### STERILISATION

Flexicare's Oxygen Bubble Tubing is supplied non-sterile.

### STORAGE CONDITIONS

Store in a cool, dry place out of direct sunlight.

### SHELF LIFE

Flexicare declares a shelf life of five years from date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product packaging.

### DISPOSAL CONSIDERATIONS

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

### PACKAGING

Primary – Polybag Polyethylene (PE)

Secondary – Carton Cardboard

#### BOX / CARTON

Product	Dimensions (mm)	Qty per carton
032-10-069	W 490mm L 320mm H 230mm	100
032-10-067	W 240mm L 190mm H 95mm	1
032-10-068	W 336mm L 310mm H 92mm	1