

## Products


[ECG](#)
[Patient Monitor](#)
[Fetal Monitor](#)
[Ultrasound](#)
[Infusion Pump & Syringe Pump](#)
[Defibrillator](#)
[Veterinary](#)
[Accessories](#)

[Online Contact](#)
[Patient Monitor](#)
[You are here](#) [Home Page](#) > [Patient Monitor](#)

### IM 12 Patient Monitor



#### Features:

- 12.1 inches high resolution color TFT display
- Optional full touch-screen for easy operation
- Attractive and durable silicon buttons with backlight
- Optional 12-lead ECG accurate analysis and further diagnostic confirmation
- S-T segment analysis and 23 types of arrhythmia analysis
- Maximum 480-hour graphic and tabular trends, 1000 NIBP records, 700 alarm events and 2-hour ECG waveform storage and review
- Large font display, OxyCRG, graphic & tabular trends, bed-bed observation and more display modes to meet different needs in clinics care
- 5 kinds of calculations: medication calculation, hemodynamic calculation, renal function calculation, oxygenation calculation and ventilation calculation
- ESU proof and defibrillation synchronization for different working environments
- Rechargeable Li-ion batteries for up to 10 hours of continuous work
- USB port enables software upgrade and also supports external USB storage
- Selectable color configurations for all parameters (32 colors)
- VGA output for external display (optional)
- Flexible mounting and trolley solution (optional)
- Wire/wireless network connect to central monitoring system (optional)

OEM Partners:



---

CopyRight ©2007-2014 Shenzhen Biocare Bio-Medical Equipment Co., Ltd. Guangdong ICP by 05081180Support: 075568.cn

Documento de muestra