

IOB WARMING UNIT

The Best Choice for Body Temperature Management to Combat Hypothermia



FDA **CE**

SAFE
RELIABLE
EASY TO USE

Rolling Stand Available

IOB MEDICAL

Copyright ©2012-2018, IOB Medical Inc. All Right Reserved.



Specifications

IOB Warming Unit



Dimensions: H28cm x D22cm x W22cm

Weight: 5.4 kg

Relative Noise Level: <55 dBA

Operating Temperature:

High: 43° ± 3°C (109.4° ± 5.4°F)

Med: 38° ± 3°C (100.4° ± 5.4°F)

Low: 32° ± 3°C (89.6° ± 5.4°F)

Device Rating:

110-120 VAC, 60 Hz, 10 A

220-240 VAC, 50 Hz, 8 A


Filter: 0.2 µm High-efficiency Filter, 500 hours

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-1-8, IEC 80601-2-35




Intraoperative Warming Blankets


Upper Body
Model-003
195cm x 64cm




XL Upper Body
Model-021
220cm x 80cm



Lower Body
Model: IOB-002
142cm x 120 cm



Torso
Model: IOB-001
142cm x 120cm



Full Body Surgical
Model: IOB-008
210cm x 120cm




Underbody Warming Blankets


Pediatric
Underbody
Model: IOB-005
91cm x 80cm




Large Pediatric
Underbody
Model: IOB-009
160cm x 80cm




Spinal
Underbody
Model: IOB-010
215cm x 100cm



Lithotomy
Underbody
Model: IOB-011
170cm x 100cm




Full Access
Underbody
Model: IOB-006
215cm x 100cm



Pediatric Warming Blankets

Pediatric
Lower Body
Model: IOB-012
142cm x 100cm



Pediatric Long
Model: IOB-014
110 x 64 cm
(200cm x 17cm)




Pediatric
Full Body
Model: IOB-007
170cm x 100cm



Pediatric
Underbody
Model: IOB-005
91cm x 80cm




Large Pediatric
Underbody
Model: IOB-009
152cm x 100cm



Postoperative Warming Blankets

Adult
Full Body
Model: IOB-004
210cm x 120cm



Pediatric
Full Body
Model: IOB-007
170cm x 100cm




Other Warming Blankets

Cath Lab
Model: IOB-015
180 x 74cm
(380cm x 17cm)



Surgical Access
Model: IOB-016
210cm x 120cm



Cardiac Access
Model: IOB-020
142cm x 120cm

